



XDR Radiology

11300 W. Olympic Blvd Ste 710, Los Angeles CA, 90064
888-XDR-XRAY (888-937-9729) Fax:(888) 937-9728
service@XDRradiology.com www.XDRradiology.com



Benefits of the XDR Anatomic Sensor

- ✓ Minimized Mesial Dead-Space
 - Improved bitewing coverage area
 - Radiographic field is moved 2mm mesial
- ✓ Beveled Corners
 - Anatomically matched to natural contours of patient's mouth
 - Broadened contact zone reduces pressure
 - More comfortable for the patient
- ✓ Smaller Sensor Size
 - The overall sensor volume is reduced by 4%
- ✓ Easier to Clean
 - No crevices to collect contaminants
- ✓ Benefits of 2 Meter Cord
 - Facilitates easier management of the cord during sensor placement
 - Reduces weight pulling on the cord for better patient comfort
 - Minimizes possible damage to the cord from feet and chair wheels
 - Reduces the chance of kinking
 - Provides for easier cleaning & storage
 - If you need to extend the length of the cord, can also provide an optional 1 meter USB extension cable at no additional cost

